

# Whole mount EdU detection of *S. mansoni*

By Jim Collins

**\*Unless otherwise noted all incubations, washes, and fixation steps should be performed on a rocker with moderate agitation.**

- 1) Label parasite with EdU (Sigma-Aldrich) either *in vitro* (10  $\mu$ M for as little as 30 minutes to 24 hours) or *in vivo* (100-200 mg/kg by IP injection into mice).
- 2) Chase EdU for desired period of time. In our experience, for *in vitro* EdU pulse periods > 4 hours can kill neoblasts, so shorted pulse periods are desired for pulse-chase experiments.
- 3) Collect mixed sexed *S. mansoni* in media (e.g., Basch Media or DMEM) with 5-10% serum. Initial steps are performed in 15ml conical tubes with 5-10ml of each solution for >50 worms or in 1.7 ml tube for <50 worms.
- 4) Separate male and female parasites by incubation in a 0.25% solution of the anesthetic ethyl 3-aminobenzoate methanesulfonate (Sigma-Aldrich, St. Louis, MO) dissolved in DMEM+FBS. Alternatively, add 1/10 volume of a 2.5% solution. Rock samples by hand gently for 1-2 minutes or until parasites are relaxed and separated.
- 5) Kill the parasites in 1ml of 0.6 M  $MgCl_2$  for ~1 min
- 6) Replace  $MgCl_2$  with 4% Formaldehyde in PBSTx, incubate 4 hours at RT
- 7) Rinse 1X with PBSTx.
- 8) Dehydrate in Methanol and store at -20°C. Samples can be stored for weeks, if not months or years, at -20°C.
- 9) Rehydrate samples in 50% Methanol solution in PBSTx, 5-10 minutes, RT
- 10) Incubate in PBSTx, 5-10 minutes, RT
- 11) Add bleaching solution (9ml  $H_2O$ , 500 $\mu$ L Formamide, 250 $\mu$ L 20x SSC, 400 $\mu$ L 30%  $H_2O_2$ ), incubate 1hr at RT under bright light.
- 12) Rinse 2x in PBSTx and then incubate in 5ug/ml Proteinase K (Invitrogen) in 1x PBSTx for 45 minutes at RT. **Note: ProK potency appears to vary greatly depending on the source and age of the enzyme. Thus, we suggest empirically determining appropriate enzyme concentration.**
- 13) Post-fix in 10 ml 4% Formaldehyde in PBSTx, 10 min at RT.
- 14) Rinse 2x with PBSTx.
- 15) Add 1ml EdU Detection solution:
  - 789  $\mu$ L 1x PBS
  - 10  $\mu$ L 100mM  $CuSO_4$
  - 1  $\mu$ L 10 mM Fluor-conjugated Azide (e.g. Azide-fluor 545, Sigma Aldrich)
  - 200  $\mu$ L 500mM Asorbic Acid (make fresh)
- 16) Incubate 30m at RT in the Dark
- 17) Wash 2-3X with PBSTx
- 18) Incubate in DAPI (1ug/ml In PBSTx) ~1hr.
- 19) Clear samples in 80% Glycerol, mount in Vectashield.

**Note:** If desired, following proteinase K incubation and post-fixation, *in situ* hybridization can be performed and EdU detection can be performed following *in situ* signal detection.